













# AUTO TRANSDUCER CABLE

<b>HEWLETT PACKARD</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0101	3	78205A,(7825 UP TO 78834) M1006A,M1006B,M1165A,M116	12		
AC0102	3	6A M1175A,M1350B,M1205A(Omnicare24C)	12		
AC0103	0,8	M1350A,M1353A,8030A,8040A Merlin,685 Viridia 24, CMS 24	12		



## Andromeda

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0201	3	Ellipse 4  Lemosa PAG MO 308 GLAC65G	8		



## Arrow

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0301	3	IABP ACAT.1 Kontron Stecker	12		
AC0302	3	Kwat II Plus  Nikolay 072454	12		



## ACTA GRAPH

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0401	3	Polygraph  Lemosa PAG MO 308 GLAC65G	7		



## BARD

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0501	3	EPU Appliance 1344;LAB.System ECG / PR.200847(Insertion)  AMP 206434.1	8		



## BBRAUN



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0601	1,5	VISICON IBP supply for CombiCable CC M010;M011;M111;M211  Lohmeier			

## BCI



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0701	3	6200;9100;9200 Advisor  Amphenol MS 3106 a.14s6.p	6		





# AUTO TRANSDUCER CABLE



<b>BMT</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0801	3	BMT.Zwonitz Universal.Adapter  RFT.BMT Connector	10		



<b>BIOMEDICAL IVY</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC0901	3	401;402;403;403D;404;405D; 700;701;704;706;707  Amphenol MS 3106 a.14s6.p	6		



<b>BIOPAC</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1001	3	DA100A; TCI100;TCI101	5		

<b>BIOTRONIK</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1101	3	EPM;HBV(20);HDR(20); Franke.Druck;MS4010;MS6010; MS8010  universal connector	5		

<b>BRUKER</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1201	0,8	Maglife MR.Compatible(another circuit)  Nicolay 072454	12		
AC1202	3	SM 784;785;786; TM910 Defigard 2002; Physiogard MB 910  Nicolay 072454	12		


<b>CARDIO CAP</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1301	3	See Datex Engstrom  Nicolay 201070	10		

<b>COBE</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1401	3	CDX III  Amphenol 3300501	4		



<b>COLIN CORPORATION</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1501	3	BP.408;CBM.3000  ITT Cannon WK6.21C	6		

# AUTO TRANSDUCER CABLE





## COROMETRIX

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1601	3	101; 110; 111; 111A; 112; 113; 115; 116; 320; 504; 506; 512; 515; 515A; 556; 556AAX; 5567AAX	12		



## CRITIKON

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1701	3	DINAMAP PLUS(8720; 8721; 8722) DINAMAP MPS Select / Portable  Amphenol MS 3106 a.14s6.p	6		









## CRITICARE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1801	3	Scholar 1 Modell 2200  Lemos PAG.MO.5Glac65g	5		
AC1802	3	Monitor 508 Monitor 1100  Amphenol MS 3106 a.14s6.p	6		



## CWE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC1901	3	PA.931;PM.1000  ITC Cannon WK6.21C	6		



## DATASCOPE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2001	3	System 80; 82; 83; 84; 90; 90T; 95; 97IABP; 850; 870; P2; P3; 2000; 2000A; 2001; 2001A; 2002; 2002A; 22001; 3000; Passport  Amphenol 3400.013	6		
AC2002	3	Expert(see Fukuda Denshi)  Fukuda Denshi	12		
AC2003	3	System built till 1980  Amphenol 3400.013	6		
AC2004	3	Point of View  Lemos PAG.MO 307 GLAC65G	7		

## DATEX



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2101	3	Cardio.Cap CC.104; CCI.104; CMO.104; AS/3; CS/3; Light Cardicap II CG.1G; CG.1GS; CG.2CS;CG.2G; CG.2GS; CH.1R; CH.1RS; CH.2; CH.2s; M.ESTP;M.ETP;M.P	10		

## DANTEC



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2201	3	21C15; 23M02; 23M04; UD5000 MK2; UD5500 MK2; Urodyne 5000; 5500; Menuet; Menuet Compact. Etcde; Duet; Video	5		

# AUTO TRANSDUCER CABLE



## DANICA

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2301	3	Dialogua 2000 Nicolay 201070 rot	10		





## DIEFENBACH

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2401	3	EM 41 X Diefenbach Connector	7		



## DISA

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2501	3	2100; 21C15 Amphenol MS 3106 a.14s6p	6		



## DRAEGER

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2601	3	Cato; Cicero EM; Vitara Parameterbox PB 8800; PM 8010; PM8014; PM 8040; PM8060; Sulla; UM3; UM3.1;	10		
AC2602	3	North American Draeger Monitor Vitalert 2000 Dash 2000 / 3000	11		





## E for M

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2701	3	see Marquette PDS 300	11		



## FUKUDA DENSHI

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2801	3	DynaScope DS 504; 831; 833; 834; 1031A; 033A; 034A; 1060; 1100; 3300; 3340; 5100; 5300; HB60; 300; 301 Fukuda Dens	12		

## GOULD

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC2901	3	13.4615, 50.52 Gould(Brush) Gould 11.5407.50	12		
AC2902	3	IM 1000 series; SM7084; 1076.1; SP 1400 itt Cannon WKS.22C	5		

## GAELTEC

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3001	3	S7B ITT Cannon WK6.21C	6		



# AUTO TRANSDUCER CABLE



<b>HELLIGE</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3101	3	Vicom(5. Generation) SMU 610; 611; 612; 613; 639; SMK 230; 231; 260; SMU,EVO Caution: Marquette Eagle	10		
AC3102	3	Servomed(4 Generation) SMV 104; 107; 111; SMS 104; 108; 109; 164; 168; 302; 304; 308 SMK 155 Hellige SMV 104	10		
AC3103	3	Midas Hellige Midas	18		
AC3104	4	Vicom Hellige SMV 104	10		
AC3105	5	Vicom Hellige SMV 104	10		
AC3106	3	3. Generation 2. Generation Amphenol 3400.013	6		
AC3107	3	Vorverstarker DC.Druck (Preamplifier DC.Pressure) Amphenol 3400.013	6		
AC3108	3	Vorverstarker DC.Druck (Preamplifier DC.Pressure) Amphenol 3400.013	6		
AC3109	3	Program 19 Amphenol 3400.013	6		
AC3110	3	Vorverstarker TF 19 (Preamplifier TF 19) Amphenol 3400.013	6		



## HONEYWELL



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3201	3	CM / RM 120.300(4,25) Honywell CM.RM 120.300	10		
AC3202	3	VR12 / VR16 / AR6 (IM; OM Series) ITT Cannon WK6.21C	6		
AC3203	3	MEDDARS 100 / 200 / 300 Series Amphenol 6706 C14.9P	9		
AC3204	3	CM 120; 130; 150; RM 300; 301 Honywell CM.RM 120.300	10		
AC3205	3	Philips Physio.Cardiopan Lemoso FGG.OB.304CLAD62	4		



# AUTO TRANSDUCER CABLE





<b>HONEYWELL</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3206	3	Philips Physio.Cardiopan Connector for Physio.Cardiopan	9		





<b>HUNGO SACHS</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3301	3	566;570 Amphenol 3400.013	6		



<b>INVIVO RESEARCH</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3401	3	3000;3200;Omni.trek TVS Millennia; Omeda; 1400 Series Amphenol MS 3106 a.14s6.p	6		

<b>Jostra</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3501	3	HLM 20 / 21 Lemos OB	5		

<b>KONE</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3601	3	M560A; M565A; M575; MA544 Nicolay 201070 rot	10		

<b>KONTRON</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3701	3	Minimon 7131 up to 7136: Fetalmon 2A; 2B Fetalogic 4000; Fetasonde 5; 5A; 5B; K2000 Balloon Pump; KAAT II; Module 10 up to 7286; Supermon Nicolay 201070 rot	12		
AC3702	3	Series 108; 200 Amphenol 3400.013	6		

<b>LOHMEIER/BBRAUN</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3801	1,5	VISICON IBP supply for ComviCable CC M010; M011;M111;M211 Lohmeier	12		
AC3802	3	M6 08; M8 08; M910; PM315 Lohmeier	12		

<b>LITTON</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC3901	3	see Hellige Servomed Hellige SMV 104	10		

# AUTO TRANSDUCER CABLE

## MARQUETTE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4001	3	2000 Series 7000 Series 7010 Series(not traded anymore) AMP 206434.1	8		
AC4002		Cath.Lab; MAC.Lab(general term) Tram Module 100 up to 600; Muse Solar 7000; 8000; 9000; 808; Cardio.Window Eagle 1000; 3000; 3100; 4000;	11		

## MEDSET

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4101	3	Mida Coronet	10		

## MENNEN

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4201	3	Horizon 1100; 2000; Series; 3000 9000 Series; XL Series Amphenol MS 3106A.14s.6p	6		
AC4202	3	EMS; ENVOY; MERCURY; SOLO; MICRO.SOLO; MR 1200 UP to 9200; Cathlab Amphenol MS 3106a.14s.6p	6		

## MESCHER

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4301	3	Yype DCB 4 B Amphenol MS 3400.013	6		

## METEK

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4401	3	DAVID Metek.Connector	15		

## MICROMEDICS

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4501	3	System 7 Urodynamic Modul66 Lemosa FGG 3P310PLWD75R	6		

## M.M.S.

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4601	3	UD 2000 Lemosa FGG.OB.304CLAD62	4		

## NDE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4701	3	Monitor Amphenol MS 3106A.14s.6p	6		

# AUTO TRANSDUCER CABLE

<b>NEC</b>					
Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4801	3	1257; 1432X; 1434X; 1439X; 2238: BioView F; Polygraph 360 Series Monitor.Connector NEC	7		

## NELLCOR PURITAN BENNET

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC4901	3	Escort 300 A Amphenol MS 3106A.14S.6P	6		
AC4902	3	PD 240; PB 245; PB245 P Nicolay 201070 rot	10		

## NEO MEDDIX

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5001	3	Urotrace Lemosa FGG.OB.304CLAD62	4		

## NIHON KOHDEN

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5101	3	LifeScope 6; 9; 12; 14; AP.311P up to 800PA; BSM.2100 up to 8500A; OEC.3200 up to 8120A; OMP.3221E up to 7203; Mercury 9510K	5		

## OLLIKONE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5201	3	533 ITT Cannon WK5 22C	5		
AC5202	3	525 ITT Cannon WK5 22C	5		

## PENTAM

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5301	3	PR22	6		

## PETAS



Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5401	3	Monitor Amphenol MS 3106A.14S.6P	6		

## PHYSIO CONTROL





Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5501	3	VSM.1; VSM.4AL VSM.5 Amphenol MS 3106A.14S.6P	6		

# AUTO TRANSDUCER CABLE



## PPG BIOMEDICAL

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5601	3	see HELLIGE VICOM Hellige SMV 104	10		



## PRUCKA

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5701	3	CLAB I (see Datascope) Amphenol 3400.013	6		
AC5702	3	CLAB II AMP 206434.1	8		



## POLYGRAPH

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5801	3	EP 12 AMP 206434.1	8		

## POLYSTAN

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC5901	3	Polystan.Monitor universal connetor	5		



## PROTOCOL

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6001	3	Propaq 104; 106; 206; 204; Encore Amphenol MS 3106A.14S.6P	6		



## PULMOKARD QUINTON

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6101	3	EpAmp; EpLap; EpXress; Q.Cath Nicolay 201070 rot	6		

## PURITAN BENNET

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6102	3	PB 224; 240; 245 Amphenol MS 3106A.14S.6P	6		



## RFT

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6301	3	BMR 1 / BMT 2 E 235; 250	10		



# AUTO TRANSDUCER CABLE

10



## ROCHE

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6401	3	220; 222; 223; 224; 7320; 7324  Amphenol 3400.013	6		















## S&W(SIMSON&WEIL)

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6501	3	System 8000 Athena; 8041; Diascope 2; Diascope Traveler; Diascope NT  Nicolay 201070 rot	10		

## SARNS/3M

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6601	3	5000; 8000; 8595; 9000 HLM  Lemosa FGG.OB.304CLAD62	5		

## SIEMENS

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6701	3	SC 6000 SC 7000	7		
AC6702	0,2	SC 6000 Adapter 7 to 10 Pin SC 7000	7		
AC6703	0,1	Sirecust 400.1281; SC 8000 Cable for Siemens Hamobox.Adapter	10		
AC6704	1,7	Sirecust 400.1281; SC 8000	10		
AC6705	5	Sirecust 400.1281; SC 8000	10		
AC6706	3	Mingograph 7(Europe); CathCor; HKL 440M; Multigraph Recor; CathLab;Micor; Cartridge; Sicor;	15		
AC6707	3	Sirecust 400.1281; SC 8000 CrosKop X-ray system; 2002 Tochtersichtgerät	10		
AC6708	0,1	Bbraun Adapter fur Siemens Hamobox 2 channels, D.Sub.Plug	15		
AC6709	3	Sirecust 300 Serie(342R)	5		
AC6710	3	Minograph 81 USA	5		

# AUTO TRANSDUCER CABLE

## SPACELABS

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6801	3	Series 400; 500; 600; 700; 1000; 3000; 90300; 90400; (90402; 90470; 90496); 90600; Alpha 7; 9; PC 1,2; UVW; 90367 portable PC Express; TEK / VIATEK; Scout	6		
AC6802	3	Module 90496 (10cm Adapter Cable)	10		
AC6803	3	Alpha 700;703 and 900	5		

## SPIEGELBERG

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC6901	3	CPP.Monitor  Lemos a FGG 3P31 OPLWD 75R	10		

## STOECKERT / SORIN

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7001	3	Stoeckert HLM CAPS Series  Amphenol T3300501	4		

## SCHUBARD

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7101	3	Schubart  Lemos a PAG MO 307 GLAC65G	5		

## SCHWARZER

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7201	3	Physioscript PD 14(since 9/98)	5		
AC7202	3	Physioscript PD 14(till 9/98) CU12.C	12		
AC7203	3	DME 47 AO; 1131A  itt Cannon WKS.22C	6		

## SCHWARZER PICKER

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7301	3	DME 48	12		

# AUTO TRANSDUCER CABLE

## SCHILLER

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7401	3	Cardioswiss CM.8; TM.7 Cardiovit AT,2S  Schiller connector	5		

## TEKTRONIX

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7501	3	413A; 414; 414.21  Amphenol MS 3106A.14S.6P	6		

## VIATEK

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7601	3	see Spacelab  Amphenol MS 3106A.14S.6P	6		

## WIEST

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7701	3	Uro.Monitor 6000.8000 Cavomat; Cavodyn  Wiet connector	12		
AC7702	3	Merkur 2000.4000  Lemosa FGG.OB.304CLAD62	7		

## KABELMANAGEMENT

Best. Nr order . NO	Lange / length	Monitor modle / Monitor model	Pins	Steckbild plug details	Steckerbild/Picture of plug
AC7801	3	Combitrans-Cable-Datascope 2-fach Combitrans-Cable-Datascope 3-fach Combitrans-Cable-Datascope 4-fach	6		
AC7802	3	Combitrans-Cable-Datex 2-fach Combitrans-Cable-Datex 3-fach Combitrans-Cable-Datex 4-fach	10		
AC7803	3	Combitrans-Cable-HP 2-fach Combitrans-Cable-HP 3-fach Combitrans-Cable-HP 4-fach	12		
AC7804	3	Combitrans-Cable-Marquette 2-fach Combitrans-Cable-Marquette 3-fach Combitrans-Cable-Marquette 4-fach	11		
AC7805	3	Combitrans-Cable-Spacelab 2-fach Combitrans-Cable-Spacelab 3-fach Combitrans-Cable-Spacelab 4-fach	6		
AC7806	3	Combitrans-Cable-Siemens 2-fach Combitrans-Cable-Siemens 3-fach Combitrans-Cable-Siemens 4-fach	10		